

OCULOPLASTIK

OPHTHALMOLOGY CATALOGUE
March 2021

Bilateral Conformers 22



NOT CE
MARKED

The most thin and flat. Sold non sterile.
19-091 [extra-small] (New single)

19-092 [small]

19-094 [regular]

Also available with 2 holes.

19-092T [small]

19-094T [regular]

Bilateral Conformers 18 and 18/20



NOT CE
MARKED

18 18/20

19-112 [s] 19-102 [s]

19-114 [reg] 19-104 [reg]

Triangular bilateral conformers with various curvatures to fit different posterior surfaces of the anophthalmic sockets. The 18/20 is a thicker and flatter model than the 18. (PMMA) Each in 2 sizes. Sold non sterile.

Sizers (Stainless steel)

With polished balls.

NOT CE
MARKED



A1-014 14mm, A1-016 16mm


A1-018 18mm, A1-019 19mm

A1-020 20mm, A1-021 21mm

A1-022 22mm

A1-SET (14,16,18,20,22)

Symblepharon Ring Conformers



NOT CE
MARKED

Ideal to avoid re-adhesion of symblepharons. (PMMA)

Round (Sold non sterile.)

21-412 [small, 20 mm]

21-414 [med., 22 mm]

21-416 [large, 24 mm]

Suction Cups



NOT CE
MARKED

For easy insertion and removal of prostheses, conformers, and ocular shields without handle. The synthetic latex-free rubber suction cups are autoclavable. 23-012LF. A dozen.

PMMA cannot be autoclaved.
Currently no devices are sold sterile.

Fine products for the
field of ophthalmology.
Other catalogues for dermatology,
plastic surgery and ophthalmology.

Oculo-Plastik, Inc.

200 Sauvé Street West, Montreal
Quebec H3L 1Y9 Canada

☎ +1 514 381-3292 • 1 888 381-3292

☎ +1 514 381-1164 • 1 800 879-1849

sales@oculoplastik.com •

www.oculoplastik.com

1p cata ofta 210303

OCULOPLASTIK

