

LASER & IPL SAFETY EYEWEAR

#31 – Small



#34 RX Insert



#33-Universal



#35- Wrap around



Most popular
#36 - Adjustable



#52 – Wrap Around Large



#39 – Large Fitover



CE CERTIFIED

	Wavelength	OD	%VLT	Color	Retail Price
Alexandrite	180-400 nm	5+			
	730-760 nm	5+	51%	Green	\$180
LIA #31,33,34,35,36,52,39	750-760 nm	7+			
Alexandrite/Diode	190-400nm	6+			
	720-810 nm	7+	33%	Pink	\$180
AXX #31,33,34,35,36,52,39	810-830 nm	5+			
Argon	190-449 nm	5+			
	450-532 nm	6+	39%	Orange	\$105
ARG #31,33,34,35,36,52,39					
Diode	190-390 nm	6+			
	785-830 nm	5+	61%	Pink	\$180
DI2 #31,33,34,35,36,52,39	800-818	6+			
HOY (Erbium)	190-400 nm	5+			
	1775-3000 nm	4+			
	1980-2550 nm	5+	41%	Grey	\$175
HOY #31,33,34,35,36,52,39	5200-11000 nm	5+			
	9000-10600 nm	6+			
YAG	190-400	5+			
	1025-1075	5+			
	1035-1070	6+	55%	Clear Grey	\$205
YG5 #31,33,34,35,36,52,39	1045-1065	7+			
Nd:YAG	190-400 nm	5+			
	1064 nm	7+			
	950-1070 nm	7+	59%	Green	\$180
YG3 #31,33,34,35,36,52,39	840-950 nm	5+			

	Wavelength	OD	%VLT	Color	Retail Price	
Alex/Diode/Yag	190-435 nm	5+				
	730-1085 nm	5+				
CYN #31,33,34,35,36,52,39	755 nm	7+	36%	Green	\$180	
	1064 nm	7+				
DYE/VBeam	576-600 nm	5+		Purple		
DYE #31,33,34,35,36,52,39	585-595 nm	6.5+	15%		\$165	
Pulsed Dye/ND Yag	190-400 nm	5+				
	582-585 nm	5+				
DYG #31,33,34,35,36,52,39	585-600 nm	6+	11%	Green	\$225	
	950-1064 nm	6+				
HYG Yag/Erbium	1025-1075 nm	5+				
	1035-1070 nm	6+				
HYG #31,33,34,35,36,52,39	1045-1065 nm	7+				
	>1775-3000 nm	4+	23%	Grey	\$210	
	1980-2550 nm	5+				
	5200-11000 nm	5+				
	9000-10600 nm	6+				
KTP/Yag/Diode	190-534 nm	7+				
	960-1064 nm	7+				
DBY #31,33,34,35,36,52,39	925-1070 nm	6+	35%	Amber	\$225	
	850-925 nm	5+				
Combination 1540/1064	190-480 nm	5+				
	800-1790 nm	5+				
IRD5 #31,33,34,35,36,52,39	820-950 nm	6+	16%	Green	\$190	
	>950-1400 nm	7+				
	>1400-1750 nm	6+				
CO2 #31,33,34,35,36,52,39	190-375 nm	5+	93%	Clear	\$55	
	10,600 nm	5+				
IPL Eyewear						
3PL #31,33,34,35,36,52,39	190-1200 nm		15%	Green	\$80	\$64

OCULOPLASTIK

LaserSecure® line of products

200 SAUVÉ West. City of MONTRÉAL, Québec, Canada H3L1Y9
 Tel (514) 381-3292 • Fax (514) 381-1164 * Toll-free tel (888) 381-3292 • Fax (800) 879-1849
www.oculoplastik.com info: sales@oculoplastik.com